



CUSHIONING, POSITIONING & SUPPORTS



REFORM SUPPORT CUSHIONS



SEAT SUPPORT



FULL BACK SUPPORT



LOW BACK SUPPORT



KNEE SUPPORT



LUMBAR SUPPORT




NECK SUPPORT



REFORM SEAT SUPPORT

Reform seat cushions combine an ergonomically designed shape with very high density, charcoal infused ActiveX™ material. This gives the user a high level of support, comfort and reduced pressure points. Anti-slip base fabric for maximum stability.


The unique molded design provides maximized support over a greater area, reducing pressure.



1. Pommel ridge to reduce adduction risk
2. Molded thigh area reducing pressure
3. Coccyx cutout to reduce pressure on tailbone area
4. Extended lower back support



PRESSURE CARE



Subject 1, Female. 128lbs Subject 2, Male. 185lbs

PRODUCT CODE

Description Code	Reform Seat Support ICR1
------------------	--------------------------

REFORM FULL BACK SUPPORT

Reform Full Back Support is designed to relieve back pressure and help posture. It naturally contours against the back with its unique shape and material. The Reform Full Back Support can be used in any type of chair or seat with an elastic strap to secure in place.

The unique shape contours to the back with extended upper back arms and two way curve for lower back support.



1. Ergonomic spine channel
2. Extended upper back and blade cushioning.
3. Rounded lumbar support
4. Targeted lower back supports



PRODUCT CODE

Description Code	Reform Full Back Support ICR2
------------------	-------------------------------

REFORM LOW BACK SUPPORT

Reform Low Back Support is designed to relieve back pressure and help posture. It naturally contours against the back with its unique shape and material. The Reform Low Back Support can be used in any type of chair or seat with an elastic strap to secure in place.

With the unique knuckle bone 2 way curve for the lower back and wider base section this support is designed to reduce lower back pressure.



1. Ergonomic spine channel
2. Rounded mid back support
3. Targeted lower back supports




PRODUCT CODE

Description	Reform Low Back Support
Code	ICR3

REFORM KNEE SUPPORT

The Reform Knee Support is designed to raise the knees or heels when in supine position. It can help to improve circulation, alleviate pain and reduce swelling or relieve pressure from bony prominences.

The unique molded design provides maximized support over a greater area reducing pressure.



1. Anti adduction ridge
2. Mounded sides securing legs in place
3. Two angled profiles for pain management and pressure relief.




PRODUCT CODE

Description	Reform Knee Support
Code	ICR4


REFORM LUMBAR

The Reform Lumbar Support can be used in Supine and Fowlers position when lying or seated. This helps to provide more lumbar and lower back support when lying down and also helps lateral positioning.

The Bed Lumbar is designed to cradle lower back, help to alleviate back pain and provide lateral positioning.



1. Thin mounded lumbar support
2. Two way angle and side support ridge



PRODUCT CODE	
Description Code	Reform Bed Lumbar ICR5

REFORM NECK SUPPORT

The Reform Neck Support helps to support the neck with unique molded shape with lateral support.

Helps to support the neck in sitting positions.




















1. Molded neck indent
2. Angled ends for lateral support



PRODUCT CODE	
Description Code	Reform Neck Support ICR6

SUPPORTA CUSHION RANGE

<p>2" Seat Cushions</p>  <p>15.75" x 15.75" 17.75" x 17.75" 17.75" x 19.75"</p>	<p>4" Seat Cushions</p>  <p>15.75" x 15.75" 17.75" x 17.75" 17.75" x 19.75"</p>	<p>Posture Wedge Cushion</p>  <p>Assists with spinal alignment</p>
<p>Bariatric Seat Cushion</p>  <p>Three layer firm cushion</p>	<p>Donut Cushion</p>  <p>For comfort while seated</p>	<p>Coccyx Cushion</p>  <p>For coccyx pain relief</p>
<p>Back Cushion</p>  <p>Relieve stiffness and pain</p>	<p>Half Lumbar Cushion</p>  <p>Supports lower back pressure</p>	<p>Knee Cushion</p>  <p>Supports pelvis and spine</p>
<p>Full Lumbar Cushion</p>  <p>Supports lower back pressure</p>	<p>Large Wedge Cushion</p>  <p>Folding cushion 26.75" x 30.75"</p>	<p>Small Wedge Cushion</p>  <p>Folding cushion 27" x 23.25"</p>
<p>Hip Abductor</p>  <p>Supports the hip - post surgery</p>	<p>Leg Raiser</p>  <p>Assists with fluid management</p>	<p>Rotator Cuff Cushion</p>  <p>Supports the shoulder - post surgery</p>
<p>Body Positioning Pillow</p>  <p>Multi-use positioning cushion</p>	<p>Body Pillow</p>  <p>Support and positioning</p>	<p>Pressure Pad</p>  <p>Multi purpose pressure relief pad</p>



VENTILATED FABRIC
Vant fabric is a three-dimensional polypropylene material that has high air permeability and is antibacterial. Using this material in the cushion covers greatly increases air movement.



PRESSURE & SUPPORT
The Icare cushions provide a huge range of solutions around positioning, alignment and pressure care.




FOR EVERYONE
Icare Supporta cushions are designed to be used at home, travelling, in the office or wherever additional comfort and support is required. Ideal for pre or post surgery relief.

SEAT CUSHION


The Icare seat cushioning range offers a range of sizes in 2in and 4in thickness. With washable zip off covers and a high density core. The Icare cushions are perfect for added comfort in wheelchairs, travelling and seating in the home.

The Supporta seat cushion range features an Visco elastic core which respond to pressure in a different way. Instead of pushing back like standard foam cushions, the Icare cushions will mold around the shape of your body taking away points of pressure.



COMMON USAGE

The Icare Seat cushions offer relief for posture-related problems, such as back pain, sciatica, and pressure points. Additionally, these cushions can aid in achieving postural alignment and pelvic stability.



ZEROTEC COVER AVAILABLE

Zerotec fabric is a medical barrier fabric that has high elasticity in every direction which reduces shear on pressure redistribution surfaces. It is waterproof and wipeable but has microporous breathable properties. This cover has welded seams and waterfall flaps for fluid control.



PRODUCT CODE


Description Code	Icare 2in Seat Cushion Standard IC540	15.75" (L) x 15.75" (W) x 2" (H)
Description Code	Icare 2in Seat Cushion Large IC545	17.75" (L) x 17.75" (W) x 2" (H)
Description Code	Icare 2in Seat Cushion Extended IC54550	17.75" (L) x 19.75" (W) x 2" (H)
Description Code	Icare 4in Seat Cushion Standard IC1040	15.75" (L) x 15.75" (W) x 4" (H)
Description Code	Icare 4in Seat Cushion Large IC1045	17.75" (L) x 17.75" (W) x 4" (H)
Description Code	Icare 4in Seat Cushion Extended IC104550	17.75" (L) x 19.75" (W) x 4" (H)



POSTURE WEDGE


The wedge shaped Icare Posture Cushion helps to restore the spine's natural shape and helps to prevent further damage to your back by removing excessive pressure on the spine. Foam insert with washable cover, wedge shape to improve posture.

Size: 15.75" (L) x 16.5" (W) x 2.75" (H) thick end



COMMON USAGE

Individuals may find relief for conditions such as, disc herniation, sciatica, lower back pain. The Posture Wedge can also assist with postural alignment and pelvic stability.



ZEROTEC COVER AVAILABLE

Zerotec fabric is a medical barrier fabric that has high elasticity in every direction which reduces shear on pressure redistribution surfaces. It is waterproof and wipeable but has microporous breathable properties. This cover has welded seams and waterfall flaps for fluid control.

PRODUCT CODE

Description Code	Icare Posture Wedge IC168
------------------	------------------------------





BARIATRIC SEAT CUSHION

Our Bariatric Seat Cushion features a 3-layer design for exceptional comfort and support. With high-density foam, ergonomic contouring, and breathable fabric, it promotes healthy posture and enhances sitting experiences for bariatric individuals.

Size: 19.5" (L) x 17.75" (W) x 4" (H)

Safe working load: 550lbs

3 Layer design: Heavy duty core



ZEROTEC COVER AVAILABLE

Zerotec fabric is a medical barrier fabric that has high elasticity in every direction which reduces shear on pressure redistribution surfaces. It is waterproof and wipeable but has microporous breathable properties. This cover has welded seams and waterfall flaps for fluid control.

PRODUCT CODE

Description Code	Icare Bariatric Seat Cushion IC188
------------------	---------------------------------------

DONUT CUSHION

This cushion is ideal for individuals who need to sit for extended periods of time and are suffering from conditions like hemorrhoids, tailbone pain, coccyx pain, sciatica or even surgeries and pre/post-natal pain. Icare Donut Cushions are designed to reduce pain and enhance your sitting comfort while at work or home. Reduces pressure in the area of discomfort.

Size: 15.75" (L) x 15.75" (W) x 2" (H) | Center hole 4.25"



COMMON USAGE

May assist users with pain relief for coccyx pain, hemorrhoid pain, postpartum recovery, pilonidal cyst or abscess or tailbone injuries or fractures.



PRODUCT CODE

Description Code	Icare Donut Cushion IC166
------------------	------------------------------

COCCYX WEDGE

Coccyx pain (Coccydynia) can be the result of a fall or a fracture. Icare Coccyx Cushions are a way to alleviate coccyx pain reducing pressure in the area of discomfort. Wedge shaped foam insert with washable cover, specifically for coccyx pain and pressure relief.

Size: 15.75" (L) x 16.5" (W) x 2.75" (H) thick end




COMMON USAGE

May assist users with pressure relief for coccyx or Coccydynia, hemorrhoid pain, postpartum recovery, pilonidal cyst or abscess or tailbone injuries or fractures.



ZEROTEC COVER AVAILABLE

Zerotec fabric is a medical barrier fabric that has high elasticity in every direction which reduces shear on pressure redistribution surfaces. It is waterproof and wipeable but has microporous breathable properties. This cover has welded seams and waterfall flaps for fluid control.



PRODUCT CODE

Description Code	Icare Coccyx Wedge IC171
------------------	-----------------------------

BACK CUSHION

The Icare Back Cushion is ergonomically designed to provide correct posture and relieve stiffness and pain. It provides targeted support and is contoured to reduce muscle strain. Pressure relieving material reduces pressure in the area of discomfort. Elastic strap with clip. Zip off washable cover.

Size: 12.5" (L) x 12.5" (W) x 3-4" (H)



COMMON USAGE

A back cushion can provide comfort and reduce pain by supporting the natural curvature of the spine. Ideal for use at home, office or travelling.




PRODUCT CODE

Description Code	Icare Back Cushion IC167
------------------	-----------------------------

HALF LUMBAR ROLL

Icare Lumbar Rolls support the natural curvature of the spinal cord relieving fatigue, pressure, and pain. Enhances correct posture, supports natural curvature of the spine. The Icare Half Lumbar Roll is small and light enough to use at home, take to the office or while travelling. Zip off washable cover.

Size: 13.75" (L) x 5.5" (W) x 3" (H)



COMMON USAGE

The Half Lumbar Roll cushion supports correct posture, neck and back alignment when seated. This cushion can also be used in bed to support the neck while lying down.



PRODUCT CODE

Description Code	Icare Half Lumbar Cushion IC170
------------------	------------------------------------

KNEE CUSHION

Icare Knee Cushion is used to alleviate back, knee, leg, and hip pain by promoting correct spinal alignment. This can also be used as a knee spacer during pregnancy, and can be placed under the legs to support the lower back and increase circulation.

High density core with zip off washable cover.

Size: 9.75" (L) x 7.75" (W) x 6" - 2.75" (D)




COMMON USAGE

The Icare Knee Cushion may assist with pelvic alignment for side sleepers. It can also reduce pain before hip or knee replacement surgery. The Knee Cushion can also be used under the ankles to off-load heel pressure.

PRODUCT CODE

Description Code	Icare Knee Cushion IC173
------------------	-----------------------------



FULL LUMBAR ROLL

Icare Lumbar Rolls support the natural curvature of the spinal cord relieving fatigue, pressure, and pain. These cushions also support correct posture, neck and back alignment when seated. The lightweight cushions are made of high-density material that retains its shape after use. Enhances correct posture, supports natural curvature of the spine. Zip off washable cover.

Size: 11.75" (L) x 4" (W) x 4" (H)




COMMON USAGE

The Full Lumbar Roll cushion can be positioned in the lumbar area to reduce pressure and support the spine.

PRODUCT CODE

Description Code	Icare Full Lumbar Roll Cushion IC169
------------------	---



LARGE BED WEDGE


An Icare Bed Wedge can be used to elevate the upper body to alleviate acid reflux, respiratory problems, and back or neck pain. It can also be used to increase circulation and reduce conditions such as varicose veins in the legs. Three in one design, the Icare bed wedge can also be placed under or on top of the mattress for alleviating head or legs. Elevation angle of 17°. Zip off washable cover.

Size: 26.75" (L) x 30.75" (W) x 8.25" (H) thick end



COMMON USAGE

The Bed Wedge can be used to elevate the head and shoulders, or the feet and lower legs.




PRODUCT CODE

Description Code	Icare Large Bed Wedge IC172
------------------	--------------------------------

SMALL BED WEDGE

An Icare Small Bed Wedge can be used to elevate the upper body to alleviate acid reflux, respiratory problems, and back or neck pain. It can also be used to increase circulation and reduce conditions such as varicose veins in the legs. Zip off washable cover. Elevation angle of 15°.

Size: 27" (L) x 23.25" (W) x 6.25" (H)



COMMON USAGE

The Bed Wedge can be used to elevate the head and shoulders, or the feet and lower legs. The smaller size is better suited to petite persons.



PRODUCT CODE

Description Code	Icare Small Bed Wedge IC182
------------------	--------------------------------

HIP ABDUCTOR CUSHION


The Hip Abductor cushion is a device used to prevent the hips from moving out of the joint. The pillow is placed between the thighs and attached to your legs with straps. Used to keep the hips stable while in bed. Adjustable velcro straps. Zip off washable cover.

Size: 17" (L) x 13" - 5" (W) x 6" (H)



COMMON USAGE

The Hip Abductor cushion may be used for immobilization or post-operative positioning and to enable healing. Also can prevent hip dislocation.



PRODUCT CODE

Description Code	Icare Hip Abductor Support IC189
------------------	----------------------------------

LEG RAISER

Icare Leg Raiser Cushion is designed to elevate the lower leg while in the supine position. With the leg cradle cut out, this assists with fluid management, positioning and pressure redistribution and can help with lower back pain.

Size: 31.5" (L) x 13.25" (W) x 9.5" (H)



COMMON USAGE

Suitable for post surgery recovery, varicose veins, swollen ankles, phlebitis oedema.



PRODUCT CODE

Description Code	Icare Leg Raiser IC184
------------------	------------------------

ROTATOR CUFF RECOVERY PILLOW


Experience optimal support and comfort for rotator cuff recovery with this specialized pillow. Its unique wedge-shaped design reduces strain, promotes proper alignment, and accelerates healing.

Size: 26.75" (L) x 30.75" (W) x 6" (H)



COMMON USAGE

The design of this support reduces weight on shoulder area for side sleepers. Used specifically for shoulder procedure recovery.




PRODUCT CODE

Description Code	Rotator Cuff Recovery Pillow IC190
------------------	---------------------------------------

BODY POSITIONING PILLOW


Designed for comfort and support during pregnancy and for positioning with our body positioning pillow. Helping with relief from pregnancy-related positioning, alleviating pressure and body alignment.

Size: 15" (L) x 9.75" (W) x 3" (H)



COMMON USAGE

This cushion can assist with positioning and pressure care for various parts of the body during sleep or rest, including spinal alignment and pregnancy support.



PRODUCT CODE

Description Code	Body Positioning Pillow IC191
------------------	----------------------------------

BODY PILLOW

The Icare Body Pillow provides all night comfort, support and positioning. With a unique hollow core fiber encasement and beaded ActiveX™ core this cushion can be molded to help a wide range of needs. Used for pregnancy, alignment, security and positioning. Helps to align spine, hips and legs to reduce pain. Zip off washable cover.

Size: 52.25" (L) x 22.5" (W) x 6.25" (H)



FIBER ENCASEMENT

The unique hollow core fiber encased cover gives the user a soft breathable comfort layer.




CHIPPED MOLDABLE CORE

ActiveX™ material in the core is chipped to provide a unique support that is able to be molded to the user.



SENSORY FRIENDLY COVER

A sensory-friendly cover features soft and gentle textures, hypoallergenic materials, and is designed to reduce noise and tactile stimulation.




PRODUCT CODE

Description Code	Icare Body Pillow ICP9
------------------	---------------------------

PRESSURE PAD

The Icare Pressure Pad can be used to provide extra pressure relief in a variety of applications when sitting or lying. The Pressure Pad can be rolled and placed under heels when sitting or lying; or provide targeted pressure relief for sacral and coccyx area. Zip off washable cover.



Size: 35.75" (L) x 18" (W) x 2" (H)



PRESSURE REDISTRIBUTION

ActiveX™ is a high grade of elastic type foam which is responsive to temperature around inflamed areas. ActiveX™ will become softer as it becomes warmer. So around a pressure point it will always soften and allow the blood to circulate freely.





NON-SLIP BASE FABRIC

The non-slip fabric on the base of the Icare Pressure Pad prevents the pad slipping when used on vinyl or slippery surfaces.



SENSORY FRIENDLY COVER

A sensory-friendly cover features soft and gentle textures, hypoallergenic materials, and is designed to reduce noise and tactile stimulation.



PRODUCT CODE

Description Code	Icare Pressure Pad IC194
------------------	-----------------------------



ICARE MEDICAL GROUP

icaremedicalgroup.com

us.support@icaremedicalgroup.com

+1 308 351 3000